

Windows® 7 (32 bit and 64 bit) Compatibility Statement

PEB XLDent is dedicated to developing the highest quality products and support for users working in current Business Editions of Microsoft operating systems, including the Windows 7™ Professional Edition, both 32 bit and 64 bit editions.

All XLDent™ Suite products are compatible with Windows 7 Professional Edition (32 bit and 64 bit). We do not recommend use, nor do we support Windows Home Premium or Home Basic editions. We strongly recommend adhering to the PEB XLDent Standard Hardware and Network Configuration Guidelines.

IMPORTANT Windows 7 (32 bit and 64 bit) Information

Be aware that not all hardware drivers are Windows 7 compatible. In addition, not all hardware drivers are 64 bit compatible. Examples of devices that may not be Windows 7 and/or 64 bit ready include printers, scanners, Digital X-ray equipment, and Intra-oral Camera units. It is the client's responsibility to verify that all software and hardware devices being utilized with XLDent Suite Products are Windows 7 compatible. In addition, it is the client's responsibility to verify that all software and hardware devices being utilized with XLDent Suite Products are 64 bit compatible in the event a Windows 7 64 bit O/S is installed. Under no circumstance will PEB XLDent be liable, nor responsible for 3rd party software applications and/or hardware components that are not Windows 7 and/or 64 bit compliant.